


How to locate your ID NOW Serial Number

The ID NOW instrument Serial number is listed in the Preferences Menu under About

Home>Preferences



The screenshot shows a list of settings: Language (English), Sound (Enabled), Brightness (50%), and About. The 'About' option is highlighted with a red box. Below the list are navigation icons: a downward arrow, a home icon, and a rightward arrow.

Preferences allow the modification of the instrument settings.


Settings include:

- Language
- Sound
- Brightness

The About screen displays information about the instrument.

Chosen options are displayed to the right of each setting.

Home>Preferences>About



The screenshot shows the 'About this Device' screen. It features the Abbott logo, the serial number 'Serial No. CB8EC413' (highlighted with a red box), and copyright information: 'Software © 2013-2019 Azzin, used under license. 2018-2019 Abbott. All rights reserved. All trademarks referenced are trademarks of either the Abbott group of companies or their respective owners.' Navigation icons at the bottom include 'Back', a home icon, and 'Info'.

Information about the instrument is displayed here.

Press Info to display further instrument information.

It is also shown on the back of your ID NOW instrument on the sticker in the bottom right corner next to this icon:



Serial Number

If you need additional information, please contact Technical Support:

Technical Support

Advice Line

Further information can be obtained from your distributor, or by contacting Technical Support on:

USA

+1 855 731 2288

ts.scr@abbott.com